TenvisHR

Activity Tracker & Heart Rate Monitor

User Manual

CAUTION

This product is not a toy. Do not allow children or pets to play with your Oaxis product. The product contains small components that can be a choking hazard.

CONSULT YOUR DOCTOR PRIOR TO USE IF YOU:

- have a medical or heart condition;
- are taking any photosensitive medicine;
- have epilepsy or are sensitive to flashing lights;
- have reduced circulation or bruise easily;
- have tendonitis, carpal tunnel syndrome, or other musculoskeletal disorders.

ACCURACY

Oaxis activity trackers are intended to be tools to provide you with information to encourage an active and healthy lifestyle. Oaxis activity trackers rely on sensors that track your movement and other metrics. The data and information provided by these devices is intended to be a close estimation of your activity and metrics tracked, but may not be completely accurate, including steps, sleep, distance, heart rate and calorie data.

WATER RESISTANT

The product is water resistant. This product is not suitable for swimming, snorkeling, diving or other water sports. Higher temperatures in places like saunas, steam rooms and hot tubs could damage the device.


APPEARANCE

1. Touch Button
2. Heart Rate Sensor
3. USB Charging Port

CHARGING THE DEVICE

1. Remove the bottom strap located at position (3).
2. Insert USB into a USB port to charge.

Note:
The strap is designed to be sturdy. To remove, carefully wriggle the band left and right.

INSTALL Tenvis HR APP

Search "Tenvis HR" in iTunes App Store or Google Play to download and install the app.

TURNING ON THE DEVICE

Press and hold the Touch Button for 3 seconds.

TURNING OFF THE DEVICE

Press the Touch Button and navigate to the 'OFF' screen. Press and Hold the Touch Button for 4 seconds. Select your option in Red and long press the Touch Button for 4 seconds to select.

DATA SYNC

Data is synced automatically when you open the Tenvis HR App. This include your phone's date and time. To sync manually, pull down the screen when you are in Tenvis HR App.

PAIRING Tenvis HR TO YOUR PHONE

1. Switch on Bluetooth in your phone's settings;
2. Allow App permissions for everything, including location.

Android: Settings > Apps & notifications > Tenvis HR > Permissions > Enable all
iOS: Health > Sources > Tenvis HR > Turn All Categories On

3. Open Tenvis HR app > click on 'Connect Tenvis HR' > Search > Select Tenvis HR

WEARING THE DEVICE


Keep A Distance
Since blood increases further up your arm, wear your tracker higher on your wrist to improve the heart rate signal during workouts. Keep a distance of about 2 to 3 finger width above your wrist bone.

Secure It
Secure the tracker to prevent it from moving. The less your tracker moves during exercise, the more accurate the heart rate tracking throughout the activity. Do not wear the band too tightly as this can restrict blood flow.

SWITCHING DISPLAY SCREEN

Tap Touch Button once to switch to navigate between different functions.

CUSTOMIZING CLOCK INTERFACE

Navigate to the 'Layouts' screen and long press the Touch Button for 4 seconds. Scroll to your desired layout and long press the Touch Button for 4 seconds to select.

CUSTOMIZING BACKGROUND IMAGE

Navigate to the Wallpaper screen and long press the Touch Button for 4 seconds. Scroll to your desired background and long press the Touch Button for 4 seconds to select.

STOPWATCH

Navigate to the 'Stopwatch' screen. Press and hold the Touch Button for 4 seconds to enter > Press the Touch Button to start and stop. To exit, press and hold Touch Button for 4 seconds.

FIND MY PHONE

Navigate to the 'Find' screen. Press and hold the Touch Button for 4 seconds to locate your smartphone. Your phone will start ringing. To exit, press and hold Touch Button for 4 seconds.

Note: Bluetooth connection must be established for this function to work.

UNPAIRING Tenvis HR FROM THE PHONE

To pair the Tenvis HR with another phone, you should first unpair it from the current phone.

Open Tenvis HR App > Me > Device Management > Connection Settings > Disconnect Tenvis HR

Note: Connect Tenvis HR band to your

smartphone before pressing on Disconnect Band to unpair it permanently.

RESETTING Tenvis HR

Open Tenvis HR App > Me > Setting > Clear Bracelet Data & App Data & Reset Bracelet

FAQ

1. **How to pair with another phone.** Before pairing with a new phone, the device has to be unpaired from its current phone.
2. **Sleep data is not recorded.** Sleep analysis is turned on at specified timing only. You may change your sleeping time in the Tenvis HR App > Me > Personal Information > Sleep Time & Wake Up Time.
3. **Calls/message notifications does not work.** Allow all permissions for the

App and allow App to work in background.

4. Heart rate value is not accurate. Secure the wristband to fit well.

Support

Visit us at <https://support.oaxis.com>

REGULATORY CONFORMANCE

Hereby, Oaxis Asia Pte Ltd declares that the radio equipment type HB3101 is in compliance with Directive 2014/53/EU.

Note: Observe the national local regulations in the location where the device is to be used. This device may be restricted for use in some or all member states of the European Union (EU).

FCC

This equipment complies with FCC RF radiation exposure limits set forth for an uncontrolled environment.

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: The manufacturer is not responsible for any radio or TV interference caused by unauthorized modifications to this equipment. Such modifications could void the user's authority to operate the equipment.

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