Infrared Digital Forehead Thermometer

User Manual

HOW TO USE

Aim towards the forehead from a distance of 5cm, press the measuring key. The temperature will display immediately. Make sure there is no hair, perspiration, cosmetic or cap covered on the forehead.

Environment Temperature

Press the measuring key, and then press LEFT button. The environment temperature will display immediately.

Memory

After measurement, press the right button to save the measurement result. Only body temperature and surface temperature can be saved.

SETTINGS

Press and hold SET button to enter the setting mode. Press SET once to switch setting from F-1 to F-5.

F-1 Units

Press LEFT to display °F. Press RIGHT to display °C.

F-2 Temperature Alarm

Press LEFT and RIGHT button to set an alarm temperature. When the body temperature exceeds a certain temperature, an alarm will be given.

Surface
Body

F-3 Modes

Press LEFT for surface temperature. Press RIGHT for body temperature.

F-4 Calibration

Default is 0. You can set the temperature deviation.

F-5: Alarm (ON-OFF)

Press LEFT button to turn on the alarm sound. Press RIGHT button to turn it off.

Measuring body temperature with forehead temperature

- Please clean your forehead before measurement to prevent hair, sweat or dirt from affecting the measurement accuracy
- Since the sensor is sensitive to the ambient temperature, do not hold the temperature sensor for a long time

Measure the temperature of objects and rooms

- After each measurement, clean the temperature head with a soft cloth and then store the machine in a dry and ventilated place
- If the product material and function break down or be damaged and affect your regular use, our company will provide a free warranty service within one year from the date of purchase.

Warranty description

- Please properly keep the warranty card as maintenance certificate;
- The warranty period is valid within one year from the date of purchase;
- During the warranty period (under the condition of normal use and maintenance), if the product failures due to the problems of its material and functions, our company will provide free maintenance and parts replacement services.

The warranty notice

For the goods beyond the warranty period or outside free maintenance service clause, our company can also provide maintenance services but need to charge the parts and maintenance fee, taking into consideration the circumstances.

The following situations are not covered by the warranty

- The equipment and components are beyond the warranty period;
- Due to human factors (such as using the product not following the instructions, error operation or improper safekeeping; human-made scratch; carrying and shock; input voltage and current not matching the product, etc.), leading to equipment failure or damage;
- Users alter or tear down the machine, or tear up the product barcode.
- OAXIS Asia Pte Ltd reserves the right of final interpretation of this warranty clause.

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